

IHE Integration Statement

Vendor	Product Name	Version	Date
Visual Medica	VM PACS	3.0.0.351	9/1/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Access to Radiology Information	Image Display	None
Access to Radiology Information	Image Manager/Archive	None
Audit Trail and Node Authentication	Secure Application	None
Consistent Time	Time Client	None
Evidence Documents	Evidence Creator	None
Evidence Documents	Image Display	None
Key Image Note	Evidence Creator	None
Key Image Note	Image Display	None
Key Image Note	Image Manager/Archive	None
Mammography Image	Image Display	None
Mammography Image	Image Manager/Archive	None
Mammography Image	Print Composer	None
NM Image	Image Display	None
NM Image	Image Manager/Archive	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Image Manager/Archive	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

HL7	http://www.hl7.org.ar/index.php/socios-institucionales
DICOM	http://visualmedica.com/vmweb/dicom.php

IHE Integration Statement

Vendor	Product Name	Version	Date
Visual Medica	VM RIS - PACS	3.0.0.351	9/1/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Scheduled Workflow	Performed Procedure Step Manager	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

HL7	http://www.hl7.org.ar/index.php/socios-institucionales
DICOM	http://visualmedica.com/vmweb/dicom.php