

IHE Integration Statement

Vendor	Product Name	Version	Date
Visual Médica	VM RIS	227	05-04-2019

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation (PIR)	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Presentation of Grouped Procedures (PGP)	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Reporting Workflow (RWF)	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Consistent Time (CT)	Time Client	None
Charge Posting (CHG)	Department System Scheduler/Order Filler	None
Audit Trail and Node Authentication (ATNA)	Secure Application	None

Internet address for vendor's IHE information: www.visualmedica.com		
Links to Standards Conformance Statements for the Implementation		
HL7: www.hl7.org DICOM: https://www.dicomstandard.org/		
Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.jira-net.or.jp/ihe-j